



- Single-use
- Designed for transesophageal transducers
- Accordion fold over applicator
- Made of soft and stretchable polyurethane material
- 12/box
- Made in the USA

The latex-free Shield® material is a proven barrier against cross-contamination and is resistant to tearing. The cover also includes an applicator which facilitates the insertion of the probe.

Item number	Description	Pckg	Price
104.0001	Sterile cover w/ sterile gel, bite block, and tape	12/box	\$213
104.0002	Non-sterile cover w/ tape	12/box	\$168
104.0003	Non-sterile cover with sterile gel, bite block, and tape	12/box	\$188